

REMOTE PATIENT MONITORING SYSTEM

Domain: Healthcare

Location: US

Project description:

Most fatal injuries happen when the patient faints/falls after getting out of bed. Insufficient lighting also causes difficulty in avoiding obstacles that could lead to falls. Even when caregivers and family members are around there could be times when the patient has to be left unattended.

- Provide comfort and peace of mind for patients having issues with falls
- Real time notification of family members and/or care providers when the patient is up
- Innovative source of lighting for patients in need
- Monitor balance, weight and patterns in movements on the mat

Requirements:

- Track patient movement when getting off the bed.
- Notify family/caregivers through alarms/app notifications if a patient fails to get back to bed within a stipulated time.
- Turn ON the lights automatically when the patient wants to get up from the bed.
- Switch off the lights automatically once the patient comes back to bed.

Solution proposed:

The core part of the project was developing:

- the switch mechanism placed within a floor mat, placed beside the bed, that could sense the patient's movement of stepping out of/in to bed.
- control device that monitors the "stepping" trigger, activates the LEDs on the mat, communicates with other "connected" devices (like lights, etc.), and sends the alerts/notifications.

The web part of the project was to provide a comprehensive interface to customize and control the various devices.

1. Profile creation
Create the customer profile. Link the product id to the customer account.
2. Manage devices:
 - a. Allows the customer to view a list of all devices that are "linked" with the controller device.
 - b. Enables the customer to configure the device settings.
 - c. Also provides a list of devices that have been "unlinked" from the controller device.
3. Enable SMS notifications by adding and verifying phone details.
4. Customize SMS notifications messages:
 - a. when the patient steps on the mat for the first time (stepping down from the bed).
 - b. when the patient steps on the mat for the second time (getting back in to bed).
 - c. as an alert if the patient fails to step on the mat the second time, which could mean the patient needs help.
5. Toggle between 24 hour SMS alerts and keeping it in "standby" for a stipulated time. This is helpful when family members/caregivers are available.
6. Advanced settings:
 - a. Set the time delay between logging in two consecutive steps on the mat. This reduces incidences of "false" triggers when the patient gets on the mat and steps on it the second time.
 - b. Specify the time delay to send the SMS after the threshold time for the patient to step back on the mat elapses.
 - c. Set the "standby" time that the mat will stay inactive after the "standby" button is pressed on the mat.
 - d. Specify the time for which the LEDs on the mat will stay ON after the patient steps on the mat for the first time.

7. Alarm settings:


- a. Set the time for which the alarm will stay active once triggered.
- b. Set the action which will trigger the alarm
- c. Set the volume level of the alarm.

Technologies used:

- Express JS framework. Node.js

Web Screens:

1. Profile




2. Linked devices



3. Controller menu





4. Controller settings



CORDLESSMAT-840D8E1C4814

SMS SETTING ADVANCED SETTING

No	Mobile No	StartTime	EndTime	First Step	Second Step	Alert	Action
1	+919820904860	00:00:00	23:59:59	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	 

Add-Modify SMS Notification :

Add & Modify Mobile No

SMS text on First Step SMS text on Second Step

SMS on Alert

Scheduling Always On (24 Hours)

5. Device menu



Welcome Anshu Singh


CORDLESSMAT-B40D8E1C4814

CORDLESSMAT-B40D8E1C4814


38E1D1278


- Settings
- Unlink Device
- Delete

6. Device settings



App Screens





The image displays three screenshots of a mobile application interface, likely for patient monitoring settings and profile management. The interface is divided into three main sections: SMS Setting, Advance Setting, and My Profile.

SMS Setting (Left Screenshot): This screen shows the configuration for adding a new SMS notification. The mobile number is +919820904860. The start time is 00:00:00 and the end time is 23:59:59. The first step is checked, and the second step is also checked. The alert is checked, and the action is set to a green checkmark and a red square.

Advance Setting (Middle Screenshot): This screen shows the configuration for various time intervals. The avoidance time is 10 seconds, the alert message time is 30 seconds, the standby time is 30 seconds, and the mat light time is 1 minute. A large blue button labeled "SAVE" is visible at the bottom.

My Profile (Right Screenshot): This screen shows the user's profile information. The first name is Anshu, the last name is Singh, and the customer account number is 0000001234. The pincode is India, the email address is farhan@embin.com, and the mobile number is +917977600586. The type is Individual, the date of birth is March 2, 1993, the country is India, and the time zone is not specified.